



# CHECK4® HCG (HL)

**Low sensitivity (1000 mIU/mL) hCG rapid test to confirm termination of pregnancy given under medical supervision**

Distributed by:

**Quadraturetech**



Product Code:

**3066-1000AG**

Size:

**1 Test**



## CHECK4® HCG (HL)

Product Code: 3066-1000AG

Check4® hCG (hl) is a rapid, qualitative test to check the outcome of first trimester Termination of Pregnancy (ToP) or miscarriage. In 5 to 10 minutes, this test will detect an hCG level of 1000 mIU/mL or more.



### TRUSTED

- Market leader chosen by renowned professionals
- Manufactured in the EU

### ACCURATE & RELIABLE


- Check4® hCG(hl) can detect hCG >1000 mIU/mL in over 99% of cases
- CE marked for home use

### QUICK

- Results in just 5 to 10 minutes

### COST EFFECTIVE

- Convenient and reduces follow up time





## HOW DOES IT WORK?

The Check4® hCG(hl) test specifically detects the hCG hormone in urine to check whether the hormone level is 1000 mIU/mL (determined against World Health Organisation reference) or greater.

If treatment for Termination of Pregnancy or miscarriage has been successful then the hCG level should drop below 1000 mIU/mL within 10 to 20 days. Check4® hCG(hl) will not detect hCG levels lower than 1000 mIU/mL, while classical pregnancy tests will be positive.

*Check4® hCG(hl) is not a classical pregnancy test and cannot be used for pregnancies where the initial hCG was lower than 1000mL.*



**For more information  
please contact:**

**Quadraturetech Diagnostics Ltd  
Station House, Cooksbridge  
Lewes, East Sussex  
BN8 4SW**

**T: 033 3321 2371  
F: 033 3321 2372  
E: [info@quadraturetech.co.uk](mailto:info@quadraturetech.co.uk)**

**[www.quadraturetech.co.uk](http://www.quadraturetech.co.uk)**

**Provider of Quality Medical Diagnostic & Research Products**

Distributed by:

**Quadraturetech**



Product Code:

**3066-1000AG**

Size:

**1 Test**

