



CHECK4[®] HCG (HL)

Low sensitivity (1000 mIU/mL) hCG rapid test to confirm termination of pregnancy given under medical supervision

Distributed by:

Quadraturetech



Product Code:

3066-1000AG

Size:

1 Test

Check4

CHECK4® HCG (HL)

Product Code: 3066-1000AG

Check4® hCG (hl) is a rapid, qualitative test to check the outcome of first trimester Termination of Pregnancy (ToP) or miscarriage. In 5 to 10 minutes, this test will detect an hCG level of 1000 mIU/mL or more.



TRUSTED

- Market leader chosen by renowned professionals
- Manufactured in the EU

QUICK

- Results in just 5 to 10 minutes

COST EFFECTIVE

- Convenient and reduces follow up time

ACCURATE & RELIABLE

- Check4® hCG (hl) can detect hCG >1000 mIU/mL in over 99% of cases
- CE marked for home use



QUADRATECH DIAGNOSTICS

Check4



HOW DOES IT WORK?

The Check4® hCG (hl) test specifically detects the hCG hormone in urine to check whether the hormone level is 1000 mIU/mL (determined against World Health Organisation reference) or greater.

If treatment for Termination of Pregnancy or miscarriage has been successful then the hCG level should drop below 1000 mIU/mL within 10 to 20 days. Check4® hCG (hl) will not detect hCG levels lower than 1000 mIU/mL, while classical pregnancy tests will be positive.

Check4® hCG(hl) is not a classical pregnancy test and cannot be used for pregnancies where the initial hCG was lower than 1000mL.



**For more information
please contact:**

**Quadraturetech Diagnostics Ltd
Station House, Cooksbridge
Lewes, East Sussex
BN8 4SW**

**T: 033 3321 2371
F: 033 3321 2372
E: info@quadraturetech.co.uk**

www.quadraturetech.co.uk

Provider of Quality Medical Diagnostic & Research Products

Distributed by:

Quadraturetech



Product Code:

3066-1000AG

Size:

1 Test

