


Check4



DISCOVER CHECK4®-hCG(hl)

Check4®-hCG(hl) is a CE-Marked low-sensitivity rapid self-test for detection of human Chorionic Gonadotropin (hCG) levels above 1000 mIU/mL in urine to determine the outcome of first-trimester termination of pregnancy (ToP) or miscarriage.



QUICK

Check4®-hCG(hl) delivers results in under 10 minutes and can be prescribed 10 – 20 days after miscarriage or medical ToP, providing results sooner than classical pregnancy tests.



ACCURATE

Trusted by leading healthcare providers in the UK, Europe, and Australia due to 100% diagnostic sensitivity and 98% diagnostic specificity whilst adhering to the highest quality standards.



EASY-TO-USE

Designed with simplicity in mind, Check4®-hCG(hl) features an intuitive design and easy-to-follow instructions, providing reliable and easy-to-interpret results.



CE-Marked for at home self-testing



Registered with the UK MHRA



Manufactured in France within the EU



CE-Marked in accordance with In-Vitro Diagnostics Directive (IVDD)



ISO 13485:2015 accredited manufacturer



Entered in the Australian Register of Therapeutic Goods (ARTG)

CHOOSE CHECK4®-hCG(hl)

As the first CE-Marked low-sensitivity hCG self- test, Check4®-hCG(hl) combines EU quality and 99% accuracy to provide a reliable solution for remote follow-up for medical abortion.

INNOVATIVE

The first low sensitivity urine hCG self-test to receive CE-marking, developed together with and for clinicians.

ACCURATE

99% Accuracy when testing on or after the day prescribed by a clinician/healthcare professional.

TELE-ABORTION

Enables remote follow-up for medical abortion, with high satisfaction reported by both patients and healthcare providers.

EU QUALITY

Made in France by ISO13485 and ISO9001 certified EU manufacturer.

TRUSTWORTHY

Trusted by leading and national healthcare providers in the UK, Europe, and Australia.

ASSURANCE

Extensively validated, CE-marked and registered with global regulators, including the UK MHRA and Australian TGA.



Looking for more information?


Visit our website for useful resources for patients and professionals

CHECK4®-hCG(hl)

[CHECK4.CO.UK](https://check4.co.uk)

CLINICAL EVIDENCE

Clinically proven, Check4®-hCG(hl) enables accurate remote self-assessment of medical abortion outcomes, aligning with international guidelines.



WHO RECOMMENDATION

The WHO recommends self-management of medical abortion in whole or in part as one of abortion care options at gestational ages <12 weeks. Self-assessment of abortion outcome using a low-sensitivity pregnancy test is as effective as assessment by a trained health worker.




COMPARABLE TO BLOOD HCG TESTS

hCG and β -hCG levels in urine have a similar or slightly higher concentration than those detected in serum. Urine tests with detection limits of 1000 mIU hCG/mL or higher have been shown to be negative 2 weeks after complete induced abortion.



COMPREHENSIVE VALIDATION

Check4®-hCG(hl) has analytical sensitivity of 1000 mIU/mL, traceable to the hCG WHO international standard, diagnostic sensitivity of 100%, diagnostic specificity of 98% and global correlation of 99% to CE-marked reference method.



PROVEN USABILITY

Healthcare providers have prescribed Check4®-hCG(hl) since 2014 to confirm ToP, supported by clinical studies showing the feasibility and high acceptance rates for 'self-assessment'.



of participants using Check4®-hCG(hl) test correctly comprehend their pregnancy status



of participants very likely/ likely to recommend the use of Check4®-hCG(hl) test




usability rated by participants on a visual analogue scale (100 very easy, 0 very difficult)

CHECK4®-hCG(hl)

[CHECK4.CO.UK](https://check4.co.uk)


HELPFUL RESOURCES

We believe in empowering both patients and healthcare professionals with the tools for better care. We have enhanced various resources including an upgraded product website, a video user guide, and instructions in multiple languages.



FOR PATIENTS

A dedicated 'For Patients' webpage providing an easy-to-use guide on how to perform the test, interpret the results, and understand the next steps.



FOR CLINICIANS

The 'For Professionals' webpage offers detailed product information, clinical evidence, latest publications, and other useful resources.



BOX LABELLING

Designed with clinicians to include a 'Perform on' date section reminding patients when to perform their test and a QR code for online resources.




VIDEO GUIDE

Product video, providing a step-by-step visual guide for the Check4®-hCG(hl) test.




MULTILINGUAL

Instructions for use available in multiple languages, supporting global patient access.



PUBLICATIONS

Library, including the latest peer-reviewed publications for the Check4®-hCG(hl) test.



For more information please contact your local distributor:

Quadrature Diagnostics - quadrature.co.uk - info@quadrature.co.uk