

Clinical Study Report of
ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)
No. CSR-SCOV2Ag–BT01/01

Product:

COVG2004001-T, Exp: 2022-04

COVG2004002-T, Exp: 2022-04

COVG2004003-T, Exp: 2022-04

Sample:

ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab) has been evaluated with specimens obtained from the patients. PCR is used as the reference method for ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab). Specimens were considered positive if PCR indicated a positive result.

Method:

The assay was performed according to the package insert. Specimens were rated as either positive or negative at 10 minutes. Do not interpret the results after 20 minutes.

All the experiments will be done in China.

Results were presented in table below

Results:

Lot# COVG2004001-T

Method		PCR		Total Results
		Positive	Negative	
ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)	Results			
	Positive	38	3	41
	Negative	2	360	362
Total Results		40	363	403

Relative Sensitivity: 95.0% (95%CI*:83.1%-99.4%)*

Relative Specificity: 99.2% (95%CI*:97.6%-99.8%)*

Relative accuracy: 98.8% (95%CI*:97.1%-99.6%)*

* Confidence Intervals



Lot# COVG2004002-T

Method		PCR		Total Results
	Results	Positive	Negative	
ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)	Positive	38	3	41
	Negative	2	360	362
	Total Results	40	363	403

Relative Sensitivity: 95.0% (95%CI*:83.1%-99.4%)*

Relative Specificity: 99.2% (95%CI*:97.6%-99.8%)*

Relative accuracy: 98.8%(95%CI*:97.1%-99.6%)*

* Confidence Intervals

Lot# COVG2004003-T

Method		PCR		Total Results
	Results	Positive	Negative	
ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)	Positive	38	3	41
	Negative	2	360	362
	Total Results	40	363	403

Relative Sensitivity: 95.0% (95%CI*:83.1%-99.4%)*

Relative Specificity: 99.2% (95%CI*:97.6%-99.8%)*

Relative accuracy: 98.8%(95%CI*:97.1%-99.6%)*

* Confidence Intervals

Conclusion:

ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab) has been calibrated against specimens that have been collected from individuals found to be either PCR positive/negative. The results were showed as blow:

Relative Sensitivity: 95.0% (95%CI*:83.1%-99.4%)*

Relative Specificity: 99.2% (95%CI*:97.6%-99.8%)*

Relative accuracy: 98.8%(95%CI*:97.1%-99.6%)*

* Confidence Intervals

For and on behalf of
ACCUBIOTECH CO., LTD.

Prepared by: Maggie Zhang Date: 2020.08.30

Review by: Andy Wang Date: 2020.08.30

Distributed By

T: 033 3321 2371
E: info@quadratediag.co.uk





Testing Data:

Results of the negative samples:

Attachment 1: the clinical data showed below that the clinical negative specimens were evaluated.

No.	Specimen	Gender	Age	PCR	ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)
1.	Nasopharyngeal Swab	Male	45	Negative	Negative
2.	Nasopharyngeal Swab	Male	57	Negative	Negative
3.	Nasopharyngeal Swab	Female	69	Negative	Negative
4.	Nasopharyngeal Swab	Female	64	Negative	Negative
5.	Nasopharyngeal Swab	Female	42	Negative	Negative
6.	Nasopharyngeal Swab	Male	52	Negative	Negative
7.	Nasopharyngeal Swab	Female	54	Negative	Negative
8.	Nasopharyngeal Swab	Male	49	Negative	Negative
9.	Nasopharyngeal Swab	Female	52	Negative	Negative
10.	Nasopharyngeal Swab	Male	71	Negative	Negative
11.	Nasopharyngeal Swab	Female	60	Negative	Negative
12.	Nasopharyngeal Swab	Male	69	Negative	Negative
13.	Nasopharyngeal Swab	Male	56	Negative	Negative
14.	Nasopharyngeal Swab	Female	64	Negative	Negative
15.	Nasopharyngeal Swab	Male	35	Negative	Negative
16.	Nasopharyngeal Swab	Male	65	Negative	Negative
17.	Nasopharyngeal Swab	Male	67	Negative	Negative
18.	Nasopharyngeal Swab	Male	50	Negative	Negative
19.	Nasopharyngeal Swab	Male	36	Negative	Negative
20.	Nasopharyngeal Swab	Male	43	Negative	Negative
21.	Nasopharyngeal Swab	Male	61	Negative	Negative
22.	Nasopharyngeal Swab	Female	56	Negative	Negative
23.	Nasopharyngeal Swab	Female	55	Negative	Negative
24.	Nasopharyngeal Swab	Female	37	Negative	Negative
25.	Nasopharyngeal Swab	Male	63	Negative	Negative
26.	Nasopharyngeal Swab	Male	49	Negative	Negative
27.	Nasopharyngeal Swab	Male	56	Negative	Negative
28.	Nasopharyngeal Swab	Female	61	Negative	Negative
29.	Nasopharyngeal Swab	Male	54	Negative	Negative
30.	Nasopharyngeal Swab	Female	56	Negative	Negative
31.	Nasopharyngeal Swab	Female	50	Negative	Negative
32.	Nasopharyngeal Swab	Male	52	Negative	Negative
33.	Nasopharyngeal Swab	Male	71	Negative	Negative
34.	Nasopharyngeal Swab	Female	72	Negative	Negative
35.	Nasopharyngeal Swab	Male	58	Negative	Negative
36.	Nasopharyngeal Swab	Male	71	Negative	Negative
37.	Nasopharyngeal Swab	Male	68	Negative	Negative
38.	Nasopharyngeal Swab	Male	54	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

39.	Nasopharyngeal Swab	Female	72	Negative	Negative
40.	Nasopharyngeal Swab	Female	45	Negative	Negative
41.	Nasopharyngeal Swab	Male	63	Negative	Negative
42.	Nasopharyngeal Swab	Male	65	Negative	Negative
43.	Nasopharyngeal Swab	Male	64	Negative	Negative
44.	Nasopharyngeal Swab	Female	68	Negative	Negative
45.	Nasopharyngeal Swab	Male	61	Negative	Negative
46.	Nasopharyngeal Swab	Male	40	Negative	Negative
47.	Nasopharyngeal Swab	Male	49	Negative	Negative
48.	Nasopharyngeal Swab	Male	64	Negative	Negative
49.	Nasopharyngeal Swab	Female	59	Negative	Negative
50.	Nasopharyngeal Swab	Female	49	Negative	Negative
51.	Nasopharyngeal Swab	Female	41	Negative	Negative
52.	Nasopharyngeal Swab	Female	36	Negative	Negative
53.	Nasopharyngeal Swab	Female	63	Negative	Negative
54.	Nasopharyngeal Swab	Male	67	Negative	Negative
55.	Nasopharyngeal Swab	Female	71	Negative	Negative
56.	Nasopharyngeal Swab	Male	38	Negative	Negative
57.	Nasopharyngeal Swab	Female	69	Negative	Negative
58.	Nasopharyngeal Swab	Female	50	Negative	Negative
59.	Nasopharyngeal Swab	Male	64	Negative	Negative
60.	Nasopharyngeal Swab	Female	37	Negative	Negative
61.	Nasopharyngeal Swab	Female	69	Negative	Negative
62.	Nasopharyngeal Swab	Female	56	Negative	Negative
63.	Nasopharyngeal Swab	Male	43	Negative	Negative
64.	Nasopharyngeal Swab	Female	71	Negative	Negative
65.	Nasopharyngeal Swab	Female	49	Negative	Negative
66.	Nasopharyngeal Swab	Male	37	Negative	Negative
67.	Nasopharyngeal Swab	Male	45	Negative	Negative
68.	Nasopharyngeal Swab	Male	44	Negative	Negative
69.	Nasopharyngeal Swab	Female	56	Negative	Negative
70.	Nasopharyngeal Swab	Female	43	Negative	Negative
71.	Nasopharyngeal Swab	Male	42	Negative	Negative
72.	Nasopharyngeal Swab	Male	64	Negative	Negative
73.	Nasopharyngeal Swab	Female	59	Negative	Negative
74.	Nasopharyngeal Swab	Male	56	Negative	Positive
75.	Nasopharyngeal Swab	Male	57	Negative	Negative
76.	Nasopharyngeal Swab	Male	39	Negative	Negative
77.	Nasopharyngeal Swab	Male	44	Negative	Negative
78.	Nasopharyngeal Swab	Female	35	Negative	Negative
79.	Nasopharyngeal Swab	Male	68	Negative	Negative
80.	Nasopharyngeal Swab	Male	50	Negative	Negative
81.	Nasopharyngeal Swab	Female	63	Negative	Negative
82.	Nasopharyngeal Swab	Male	37	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

83.	Nasopharyngeal Swab	Male	68	Negative	Negative
84.	Nasopharyngeal Swab	Female	59	Negative	Positive
85.	Nasopharyngeal Swab	Female	65	Negative	Negative
86.	Nasopharyngeal Swab	Male	72	Negative	Negative
87.	Nasopharyngeal Swab	Male	53	Negative	Negative
88.	Nasopharyngeal Swab	Male	45	Negative	Negative
89.	Nasopharyngeal Swab	Male	45	Negative	Negative
90.	Nasopharyngeal Swab	Female	49	Negative	Negative
91.	Nasopharyngeal Swab	Male	46	Negative	Negative
92.	Nasopharyngeal Swab	Male	35	Negative	Negative
93.	Nasopharyngeal Swab	Male	45	Negative	Negative
94.	Nasopharyngeal Swab	Female	66	Negative	Negative
95.	Nasopharyngeal Swab	Male	40	Negative	Negative
96.	Nasopharyngeal Swab	Female	68	Negative	Negative
97.	Nasopharyngeal Swab	Female	42	Negative	Negative
98.	Nasopharyngeal Swab	Male	53	Negative	Negative
99.	Nasopharyngeal Swab	Female	35	Negative	Negative
100.	Nasopharyngeal Swab	Female	67	Negative	Negative
101.	Nasopharyngeal Swab	Female	70	Negative	Negative
102.	Nasopharyngeal Swab	Male	51	Negative	Negative
103.	Nasopharyngeal Swab	Female	72	Negative	Negative
104.	Nasopharyngeal Swab	Female	58	Negative	Negative
105.	Nasopharyngeal Swab	Male	62	Negative	Negative
106.	Nasopharyngeal Swab	Female	53	Negative	Negative
107.	Nasopharyngeal Swab	Female	38	Negative	Negative
108.	Nasopharyngeal Swab	Female	41	Negative	Negative
109.	Nasopharyngeal Swab	Female	46	Negative	Negative
110.	Nasopharyngeal Swab	Female	66	Negative	Negative
111.	Nasopharyngeal Swab	Female	47	Negative	Negative
112.	Nasopharyngeal Swab	Female	67	Negative	Negative
113.	Nasopharyngeal Swab	Male	68	Negative	Negative
114.	Nasopharyngeal Swab	Male	43	Negative	Negative
115.	Nasopharyngeal Swab	Female	50	Negative	Negative
116.	Nasopharyngeal Swab	Male	38	Negative	Negative
117.	Nasopharyngeal Swab	Female	37	Negative	Negative
118.	Nasopharyngeal Swab	Male	38	Negative	Negative
119.	Nasopharyngeal Swab	Female	66	Negative	Negative
120.	Nasopharyngeal Swab	Male	35	Negative	Negative
121.	Nasopharyngeal Swab	Male	42	Negative	Negative
122.	Nasopharyngeal Swab	Male	54	Negative	Negative
123.	Nasopharyngeal Swab	Male	70	Negative	Negative
124.	Nasopharyngeal Swab	Male	57	Negative	Negative
125.	Nasopharyngeal Swab	Female	37	Negative	Negative
126.	Nasopharyngeal Swab	Female	36	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

127.	Nasopharyngeal Swab	Male	54	Negative	Negative
128.	Nasopharyngeal Swab	Male	56	Negative	Negative
129.	Nasopharyngeal Swab	Male	56	Negative	Negative
130.	Nasopharyngeal Swab	Female	71	Negative	Negative
131.	Nasopharyngeal Swab	Female	52	Negative	Negative
132.	Nasopharyngeal Swab	Male	69	Negative	Negative
133.	Nasopharyngeal Swab	Male	41	Negative	Negative
134.	Nasopharyngeal Swab	Male	36	Negative	Negative
135.	Nasopharyngeal Swab	Female	70	Negative	Negative
136.	Nasopharyngeal Swab	Male	69	Negative	Negative
137.	Nasopharyngeal Swab	Male	49	Negative	Negative
138.	Nasopharyngeal Swab	Female	68	Negative	Negative
139.	Nasopharyngeal Swab	Male	35	Negative	Negative
140.	Nasopharyngeal Swab	Female	69	Negative	Negative
141.	Nasopharyngeal Swab	Female	41	Negative	Negative
142.	Nasopharyngeal Swab	Female	49	Negative	Negative
143.	Nasopharyngeal Swab	Male	37	Negative	Negative
144.	Nasopharyngeal Swab	Female	61	Negative	Negative
145.	Nasopharyngeal Swab	Male	60	Negative	Negative
146.	Nasopharyngeal Swab	Female	56	Negative	Negative
147.	Nasopharyngeal Swab	Male	49	Negative	Negative
148.	Nasopharyngeal Swab	Female	36	Negative	Negative
149.	Nasopharyngeal Swab	Female	64	Negative	Negative
150.	Nasopharyngeal Swab	Male	58	Negative	Negative
151.	Nasopharyngeal Swab	Male	67	Negative	Negative
152.	Nasopharyngeal Swab	Male	50	Negative	Negative
153.	Nasopharyngeal Swab	Male	54	Negative	Negative
154.	Nasopharyngeal Swab	Female	35	Negative	Negative
155.	Nasopharyngeal Swab	Male	47	Negative	Negative
156.	Nasopharyngeal Swab	Male	51	Negative	Negative
157.	Nasopharyngeal Swab	Female	58	Negative	Negative
158.	Nasopharyngeal Swab	Male	41	Negative	Negative
159.	Nasopharyngeal Swab	Male	38	Negative	Negative
160.	Nasopharyngeal Swab	Male	41	Negative	Negative
161.	Nasopharyngeal Swab	Male	39	Negative	Negative
162.	Nasopharyngeal Swab	Male	44	Negative	Negative
163.	Nasopharyngeal Swab	Male	49	Negative	Negative
164.	Nasopharyngeal Swab	Male	42	Negative	Negative
165.	Nasopharyngeal Swab	Female	44	Negative	Negative
166.	Nasopharyngeal Swab	Male	35	Negative	Negative
167.	Nasopharyngeal Swab	Female	41	Negative	Negative
168.	Nasopharyngeal Swab	Female	44	Negative	Negative
169.	Nasopharyngeal Swab	Male	56	Negative	Negative
170.	Nasopharyngeal Swab	Male	55	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

171.	Nasopharyngeal Swab	Female	39	Negative	Negative
172.	Nasopharyngeal Swab	Male	49	Negative	Negative
173.	Nasopharyngeal Swab	Female	39	Negative	Negative
174.	Nasopharyngeal Swab	Female	44	Negative	Negative
175.	Nasopharyngeal Swab	Male	59	Negative	Negative
176.	Nasopharyngeal Swab	Female	56	Negative	Negative
177.	Nasopharyngeal Swab	Male	44	Negative	Negative
178.	Nasopharyngeal Swab	Male	50	Negative	Negative
179.	Nasopharyngeal Swab	Male	36	Negative	Negative
180.	Nasopharyngeal Swab	Male	47	Negative	Negative
181.	Nasopharyngeal Swab	Male	38	Negative	Negative
182.	Nasopharyngeal Swab	Female	52	Negative	Negative
183.	Nasopharyngeal Swab	Male	36	Negative	Negative
184.	Nasopharyngeal Swab	Female	43	Negative	Negative
185.	Nasopharyngeal Swab	Female	59	Negative	Negative
186.	Nasopharyngeal Swab	Male	68	Negative	Negative
187.	Nasopharyngeal Swab	Male	60	Negative	Negative
188.	Nasopharyngeal Swab	Male	57	Negative	Negative
189.	Nasopharyngeal Swab	Male	45	Negative	Negative
190.	Nasopharyngeal Swab	Female	70	Negative	Negative
191.	Nasopharyngeal Swab	Female	69	Negative	Negative
192.	Nasopharyngeal Swab	Male	72	Negative	Negative
193.	Nasopharyngeal Swab	Female	72	Negative	Negative
194.	Nasopharyngeal Swab	Female	41	Negative	Negative
195.	Nasopharyngeal Swab	Female	61	Negative	Negative
196.	Nasopharyngeal Swab	Male	54	Negative	Negative
197.	Nasopharyngeal Swab	Female	53	Negative	Negative
198.	Nasopharyngeal Swab	Female	71	Negative	Negative
199.	Nasopharyngeal Swab	Male	39	Negative	Negative
200.	Nasopharyngeal Swab	Female	65	Negative	Negative
201.	Nasopharyngeal Swab	Male	40	Negative	Negative
202.	Nasopharyngeal Swab	Female	46	Negative	Negative
203.	Nasopharyngeal Swab	Female	70	Negative	Negative
204.	Nasopharyngeal Swab	Male	52	Negative	Negative
205.	Nasopharyngeal Swab	Male	57	Negative	Negative
206.	Nasopharyngeal Swab	Male	54	Negative	Negative
207.	Nasopharyngeal Swab	Female	41	Negative	Negative
208.	Nasopharyngeal Swab	Female	56	Negative	Negative
209.	Nasopharyngeal Swab	Male	57	Negative	Negative
210.	Nasopharyngeal Swab	Male	71	Negative	Negative
211.	Nasopharyngeal Swab	Female	45	Negative	Negative
212.	Nasopharyngeal Swab	Female	57	Negative	Negative
213.	Nasopharyngeal Swab	Female	35	Negative	Negative
214.	Nasopharyngeal Swab	Female	58	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

215.	Nasopharyngeal Swab	Female	67	Negative	Negative
216.	Nasopharyngeal Swab	Male	35	Negative	Negative
217.	Nasopharyngeal Swab	Male	42	Negative	Positive
218.	Nasopharyngeal Swab	Male	36	Negative	Negative
219.	Nasopharyngeal Swab	Female	60	Negative	Negative
220.	Nasopharyngeal Swab	Female	68	Negative	Negative
221.	Nasopharyngeal Swab	Male	66	Negative	Negative
222.	Nasopharyngeal Swab	Male	49	Negative	Negative
223.	Nasopharyngeal Swab	Female	56	Negative	Negative
224.	Nasopharyngeal Swab	Female	36	Negative	Negative
225.	Nasopharyngeal Swab	Female	35	Negative	Negative
226.	Nasopharyngeal Swab	Female	72	Negative	Negative
227.	Nasopharyngeal Swab	Female	39	Negative	Negative
228.	Nasopharyngeal Swab	Female	50	Negative	Negative
229.	Nasopharyngeal Swab	Male	72	Negative	Negative
230.	Nasopharyngeal Swab	Male	64	Negative	Negative
231.	Nasopharyngeal Swab	Female	43	Negative	Negative
232.	Nasopharyngeal Swab	Male	60	Negative	Negative
233.	Nasopharyngeal Swab	Male	68	Negative	Negative
234.	Nasopharyngeal Swab	Male	67	Negative	Negative
235.	Nasopharyngeal Swab	Female	63	Negative	Negative
236.	Nasopharyngeal Swab	Female	53	Negative	Negative
237.	Nasopharyngeal Swab	Female	54	Negative	Negative
238.	Nasopharyngeal Swab	Female	53	Negative	Negative
239.	Nasopharyngeal Swab	Female	39	Negative	Negative
240.	Nasopharyngeal Swab	Male	59	Negative	Negative
241.	Nasopharyngeal Swab	Male	48	Negative	Negative
242.	Nasopharyngeal Swab	Female	44	Negative	Negative
243.	Nasopharyngeal Swab	Male	70	Negative	Negative
244.	Nasopharyngeal Swab	Female	52	Negative	Negative
245.	Nasopharyngeal Swab	Female	54	Negative	Negative
246.	Nasopharyngeal Swab	Male	67	Negative	Negative
247.	Nasopharyngeal Swab	Male	36	Negative	Negative
248.	Nasopharyngeal Swab	Male	35	Negative	Negative
249.	Nasopharyngeal Swab	Female	36	Negative	Negative
250.	Nasopharyngeal Swab	Female	57	Negative	Negative
251.	Nasopharyngeal Swab	Female	68	Negative	Negative
252.	Nasopharyngeal Swab	Female	63	Negative	Negative
253.	Nasopharyngeal Swab	Male	49	Negative	Negative
254.	Nasopharyngeal Swab	Female	60	Negative	Negative
255.	Nasopharyngeal Swab	Male	68	Negative	Negative
256.	Nasopharyngeal Swab	Female	61	Negative	Negative
257.	Nasopharyngeal Swab	Female	39	Negative	Negative
258.	Nasopharyngeal Swab	Male	43	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

259.	Nasopharyngeal Swab	Female	59	Negative	Negative
260.	Nasopharyngeal Swab	Male	66	Negative	Negative
261.	Nasopharyngeal Swab	Male	52	Negative	Negative
262.	Nasopharyngeal Swab	Female	63	Negative	Negative
263.	Nasopharyngeal Swab	Male	65	Negative	Negative
264.	Nasopharyngeal Swab	Male	52	Negative	Negative
265.	Nasopharyngeal Swab	Male	44	Negative	Negative
266.	Nasopharyngeal Swab	Female	48	Negative	Negative
267.	Nasopharyngeal Swab	Female	64	Negative	Negative
268.	Nasopharyngeal Swab	Female	52	Negative	Negative
269.	Nasopharyngeal Swab	Female	68	Negative	Negative
270.	Nasopharyngeal Swab	Male	48	Negative	Negative
271.	Nasopharyngeal Swab	Female	43	Negative	Negative
272.	Nasopharyngeal Swab	Female	69	Negative	Negative
273.	Nasopharyngeal Swab	Male	51	Negative	Negative
274.	Nasopharyngeal Swab	Female	58	Negative	Negative
275.	Nasopharyngeal Swab	Female	40	Negative	Negative
276.	Nasopharyngeal Swab	Female	63	Negative	Negative
277.	Nasopharyngeal Swab	Male	71	Negative	Negative
278.	Nasopharyngeal Swab	Male	61	Negative	Negative
279.	Nasopharyngeal Swab	Male	54	Negative	Negative
280.	Nasopharyngeal Swab	Male	70	Negative	Negative
281.	Nasopharyngeal Swab	Male	71	Negative	Negative
282.	Nasopharyngeal Swab	Male	47	Negative	Negative
283.	Nasopharyngeal Swab	Male	66	Negative	Negative
284.	Nasopharyngeal Swab	Female	36	Negative	Negative
285.	Nasopharyngeal Swab	Male	60	Negative	Negative
286.	Nasopharyngeal Swab	Male	47	Negative	Negative
287.	Nasopharyngeal Swab	Male	53	Negative	Negative
288.	Nasopharyngeal Swab	Female	56	Negative	Negative
289.	Nasopharyngeal Swab	Female	67	Negative	Negative
290.	Nasopharyngeal Swab	Male	71	Negative	Negative
291.	Nasopharyngeal Swab	Female	61	Negative	Negative
292.	Nasopharyngeal Swab	Male	71	Negative	Negative
293.	Nasopharyngeal Swab	Female	61	Negative	Negative
294.	Nasopharyngeal Swab	Female	49	Negative	Negative
295.	Nasopharyngeal Swab	Female	64	Negative	Negative
296.	Nasopharyngeal Swab	Male	55	Negative	Negative
297.	Nasopharyngeal Swab	Female	50	Negative	Negative
298.	Nasopharyngeal Swab	Female	71	Negative	Negative
299.	Nasopharyngeal Swab	Male	70	Negative	Negative
300.	Nasopharyngeal Swab	Female	53	Negative	Negative
301.	Nasopharyngeal Swab	Male	61	Negative	Negative
302.	Nasopharyngeal Swab	Female	42	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

303.	Nasopharyngeal Swab	Male	58	Negative	Negative
304.	Nasopharyngeal Swab	Female	46	Negative	Negative
305.	Nasopharyngeal Swab	Male	44	Negative	Negative
306.	Nasopharyngeal Swab	Male	41	Negative	Negative
307.	Nasopharyngeal Swab	Male	54	Negative	Negative
308.	Nasopharyngeal Swab	Female	51	Negative	Negative
309.	Nasopharyngeal Swab	Male	53	Negative	Negative
310.	Nasopharyngeal Swab	Male	54	Negative	Negative
311.	Nasopharyngeal Swab	Male	36	Negative	Negative
312.	Nasopharyngeal Swab	Female	61	Negative	Negative
313.	Nasopharyngeal Swab	Female	44	Negative	Negative
314.	Nasopharyngeal Swab	Male	52	Negative	Negative
315.	Nasopharyngeal Swab	Female	53	Negative	Negative
316.	Nasopharyngeal Swab	Female	35	Negative	Negative
317.	Nasopharyngeal Swab	Female	44	Negative	Negative
318.	Nasopharyngeal Swab	Female	56	Negative	Negative
319.	Nasopharyngeal Swab	Female	57	Negative	Negative
320.	Nasopharyngeal Swab	Female	52	Negative	Negative
321.	Nasopharyngeal Swab	Male	36	Negative	Negative
322.	Nasopharyngeal Swab	Female	69	Negative	Negative
323.	Nasopharyngeal Swab	Male	37	Negative	Negative
324.	Nasopharyngeal Swab	Male	54	Negative	Negative
325.	Nasopharyngeal Swab	Female	59	Negative	Negative
326.	Nasopharyngeal Swab	Male	70	Negative	Negative
327.	Nasopharyngeal Swab	Female	72	Negative	Negative
328.	Nasopharyngeal Swab	Male	63	Negative	Negative
329.	Nasopharyngeal Swab	Male	71	Negative	Negative
330.	Nasopharyngeal Swab	Female	61	Negative	Negative
331.	Nasopharyngeal Swab	Female	35	Negative	Negative
332.	Nasopharyngeal Swab	Male	63	Negative	Negative
333.	Nasopharyngeal Swab	Male	59	Negative	Negative
334.	Nasopharyngeal Swab	Female	57	Negative	Negative
335.	Nasopharyngeal Swab	Female	66	Negative	Negative
336.	Nasopharyngeal Swab	Female	51	Negative	Negative
337.	Nasopharyngeal Swab	Female	54	Negative	Negative
338.	Nasopharyngeal Swab	Male	64	Negative	Negative
339.	Nasopharyngeal Swab	Female	59	Negative	Negative
340.	Nasopharyngeal Swab	Female	45	Negative	Negative
341.	Nasopharyngeal Swab	Male	46	Negative	Negative
342.	Nasopharyngeal Swab	Female	57	Negative	Negative
343.	Nasopharyngeal Swab	Male	50	Negative	Negative
344.	Nasopharyngeal Swab	Female	52	Negative	Negative
345.	Nasopharyngeal Swab	Male	69	Negative	Negative
346.	Nasopharyngeal Swab	Male	66	Negative	Negative



AccuBioTech Co.,Ltd

Accurate, Reliable, Cost Effective

347.	Nasopharyngeal Swab	Female	53	Negative	Negative
348.	Nasopharyngeal Swab	Female	51	Negative	Negative
349.	Nasopharyngeal Swab	Male	71	Negative	Negative
350.	Nasopharyngeal Swab	Female	64	Negative	Negative
351.	Nasopharyngeal Swab	Female	46	Negative	Negative
352.	Nasopharyngeal Swab	Male	52	Negative	Negative
353.	Nasopharyngeal Swab	Female	55	Negative	Negative
354.	Nasopharyngeal Swab	Male	63	Negative	Negative
355.	Nasopharyngeal Swab	Female	60	Negative	Negative
356.	Nasopharyngeal Swab	Male	72	Negative	Negative
357.	Nasopharyngeal Swab	Female	37	Negative	Negative
358.	Nasopharyngeal Swab	Male	38	Negative	Negative
359.	Nasopharyngeal Swab	Male	52	Negative	Negative
360.	Nasopharyngeal Swab	Female	68	Negative	Negative
361.	Nasopharyngeal Swab	Male	67	Negative	Negative
362.	Nasopharyngeal Swab	Female	44	Negative	Negative
363.	Nasopharyngeal Swab	Female	66	Negative	Negative



Results of positive samples

Attachment 2: the clinical data showed below that the clinical positive specimens were evaluated.

No.	Specimen	Gender	Age	Days from onset of Symptoms	PCR	ACCU-TELL® SARS-CoV-2 Ag Cassette (Nasopharyngeal Swab)
1	Nasopharyngeal Swab	Female	50	1	Positive	Positive
2	Nasopharyngeal Swab	Male	42	0	Positive	Positive
3	Nasopharyngeal Swab	Female	20	1	Positive	Positive
4	Nasopharyngeal Swab	Male	38	1	Positive	Positive
5	Nasopharyngeal Swab	Female	39	4	Positive	Positive
6	Nasopharyngeal Swab	Female	34	2	Positive	Positive
7	Nasopharyngeal Swab	Female	51	2	Positive	Positive
8	Nasopharyngeal Swab	Female	67	1	Positive	Positive
9	Nasopharyngeal Swab	Male	70	1	Positive	Positive
10	Nasopharyngeal Swab	Female	27	2	Positive	Positive
11	Nasopharyngeal Swab	Male	54	1	Positive	Positive
12	Nasopharyngeal Swab	Female	22	3	Positive	Positive
13	Nasopharyngeal Swab	Female	36	5	Positive	Positive
14	Nasopharyngeal Swab	Female	47	2	Positive	Positive
15	Nasopharyngeal Swab	Male	62	5	Positive	Positive
16	Nasopharyngeal Swab	Female	52	6	Positive	Positive
17	Nasopharyngeal Swab	Male	41	2	Positive	Positive
18	Nasopharyngeal Swab	Male	49	1	Positive	Positive
19	Nasopharyngeal Swab	Female	50	3	Positive	Positive
20	Nasopharyngeal Swab	Male	63	2	Positive	Positive
21	Nasopharyngeal Swab	Male	16	6	Positive	Positive
22	Nasopharyngeal Swab	Female	69	5	Positive	Positive
23	Nasopharyngeal Swab	Female	55	1	Positive	Positive
24	Nasopharyngeal Swab	Female	34	3	Positive	Positive
25	Nasopharyngeal Swab	Female	28	0	Positive	Positive
26	Nasopharyngeal Swab	Male	36	2	Positive	Positive
27	Nasopharyngeal Swab	Female	70	1	Positive	Positive
28	Nasopharyngeal Swab	Female	62	0	Positive	Positive
29	Nasopharyngeal Swab	Female	13	3	Positive	Positive
30	Nasopharyngeal Swab	Female	69	3	Positive	Negative
31	Nasopharyngeal Swab	Male	58	2	Positive	Positive
32	Nasopharyngeal Swab	Female	35	1	Positive	Positive
33	Nasopharyngeal Swab	Female	41	2	Positive	Positive
34	Nasopharyngeal Swab	Male	38	2	Positive	Positive
35	Nasopharyngeal Swab	Female	42	1	Positive	Positive
36	Nasopharyngeal Swab	Female	53	3	Positive	Positive
37	Nasopharyngeal Swab	Male	62	2	Positive	Positive
38	Nasopharyngeal Swab	Male	29	0	Positive	Positive
39	Nasopharyngeal Swab	Male	33	1	Positive	Positive
40	Nasopharyngeal Swab	Male	46	2	Positive	Negative