



# COVID-19 RAPID NASAL SWAB ANTIGEN TEST

A rapid in vitro immunochromatographic assay  
for the qualitative detection of SARS-CoV-2  
nucleocapsid protein antigen.

Distributed by:

**Quadraturetech**

Product Code:

**M200969**

Registered With:



Federal Institute  
for Drugs  
and Medical Devices



Size:

**10 Tests**

## COVID-19 RAPID ANTIGEN TEST

Product Code: M200969

The Covid-19 Rapid Antigen Test is a quick, reliable and cost effective way of detecting active COVID-19 infections.

It is ideal for employers wanting to diagnose potential infections within their workforce.



### TRUSTED

- Supplier to the NHS for 25+ years
- EU Manufactured
- Evaluated in Germany and USA
- Registered with German Federal Institute for Drugs and Medical Devices (BfArM)

### QUICK & COST EFFECTIVE

- Results in just 15 minutes
- All equipment provided
- Competitively priced


### ACCURATE & RELIABLE

- Market leader with 98.74% accuracy vs PCR Testing in clinical trials
- CE marked for professional use
- ISO 13485 certification



[www.quadratech.co.uk](http://www.quadratech.co.uk)

## QUADRATECH DIAGNOSTICS



### WHAT IS ANTIGEN TESTING?

It is a new form of testing for Covid-19. Rapid antigen testing represents a major advancement in helping to bring the pandemic under control.

Antigen tests look for pieces of proteins that make up the SARS-CoV-2 virus, helping to determine if the person has an active Covid-19 infection.

The test will detect active Covid-19 infections in both asymptomatic and symptomatic patients.


*For medical healthcare professional use only - the physician should evaluate the data obtained using this kit in the light of the other clinical information available.*

[www.quadratech.co.uk](http://www.quadratech.co.uk)



# COVID-19 RAPID ANTIGEN TEST

Product Code: M200969



## WHY BUY ANTIGEN TESTS?

Our test will detect Covid-19 in both asymptomatic and symptomatic patients.

Performing the test is simple and does not require expensive equipment to read results.

Test results are delivered on-site in 15 mins. No need to send samples away for analysis.

Unlike other antigen tests, our test looks for nucleocapsid protein making it the gold standard for rapid antigen testing.

Highly accurate with 99.22% Specificity, 96.72% Sensitivity & 98.74% Accuracy during evaluation conducted in USA


[www.quadratech.co.uk](http://www.quadratech.co.uk)

QUADRATECH DIAGNOSTICS

## HOW DOES IT WORK?



1. A nasal swab is used to obtain a sample for testing.
2. The Swab is placed in the extraction tube to incubate for 1 minute.
3. The extracted solution is used as the sample material for the test.
4. 4 drops of the test sample are added to the testing cartridge.
5. After 15 minutes, the test strip shows a positive or negative result.



[www.quadratech.co.uk](http://www.quadratech.co.uk)



**For more information  
please contact:**

**Quadratesh Diagnostics Ltd  
Station House, Cooksbridge  
Lewes, East Sussex  
BN8 4SW**

**T: 033 3321 2371  
F: 033 3321 2372  
E: [info@quadratesh.co.uk](mailto:info@quadratesh.co.uk)**

**[www.quadratesh.co.uk](http://www.quadratesh.co.uk)**

**Provider of Quality Medical Diagnostic & Research Products**

Distributed by:

**Quadratesh**

Product Code:

**M200969**

Registered With:



**Federal Institute  
for Drugs  
and Medical Devices**



Size:

**10 Tests**