



**For more information  
please contact:**

Quadrtech Diagnostics Ltd  
P.O. Box 167  
Epsom  
Surrey KT18 7YL

P: 033 3321 2371  
F: 033 3321 2372  
E: info@quadrtech.co.uk

[www.quadrtech.co.uk](http://www.quadrtech.co.uk)

**Quality products for diagnostic & research testing**



P: 033 3321 2371  
F: 033 3321 2372  
E: info@quadrtech.co.uk

[www.quadrtech.co.uk](http://www.quadrtech.co.uk)



# CALIBRATORS & CONTROLS FOR ANTITHROMBOTIC DRUGS

For the measurement of proteins involved in regulation of  
blood haemostasis or thrombosis and anticoagulant drugs

Distributed by:  
**Quadrtech**



Product Code:  
**Various**

Size:  
**See inside for  
details**

Distributed by:  
**Quadrtech**



Product Code:  
**Various**

Size:  
**See inside for  
details**

## CALIBRATORS & CONTROLS

Product Code: Various

Quadratesch Diagnostics offers an extensive range of anticoagulant calibrators and controls, including the various forms of Heparin, hirudin-like substances and the new DOACS such as Rivaroxaban or Edoxaban.

The complete line of calibrators and controls are standardized against WHO International Standards (when available) for the measurement of the various proteins involved in regulation of blood haemostasis and thrombosis as well as for anticoagulant drugs.



### CONVENIENT & ECONOMIC SIZES

- Controls are offered in packs containing 2 levels of freeze-dried material each in 6 vials of 1ml to minimise waste
- Calibrators are usually 3 levels each in 4 vials of 1ml for ease of curve construction and also to minimise waste

### GOOD PATIENT MANAGEMENT

- To aid drug management, some of the new DOACS are available as low-level calibrators and controls, such as Edoxaban and Apixaban

### CE MARKED

- All materials are CE marked and suitable for use with chromogenic based drug assay kits on many automates

### ELIMINATING THE EFFECTS OF DOACS

- To eliminate drug interference in patient plasmas prior to testing, DOAC-Remove and DOAC-Stop tablets are now also available

## CALIBRATOR & CONTROL PRODUCTS

Code	Product	Size
222001	BIOPHEN™ Heparin Calibrator	5 x 4 x 1 ml
222201	BIOPHEN™ Orgaran® Calibrator	5 x 4 x 1 ml
222301	BIOPHEN™ UFH Calibrator	5 x 4 x 1 ml
222501	BIOPHEN™ Arixtra® Calibrator	4 x 3 x 1 ml
222701	BIOPHEN™ Rivaroxaban Plasma Calibrator	3 x 4 x 1 ml
222801	Dabigatran Plasma Calibrator	3 x 4 x 1 ml
222901	BIOPHEN™ Dabigatran Calibrator Low	3 x 4 x 1 ml
223001	BIOPHEN™ LMWH Control Plasma	2 x 6 x 1 ml
223101	BIOPHEN™ UFH Control Plasma	2 x 6 x 1 ml
223501	BIOPHEN™ Orgaran® Control	2 x 6 x 1 ml
223701	BIOPHEN™ LMWH Control Low	2 x 6 x 1 ml
223801	BIOPHEN™ LMWH Control C3	12 x 1 ml
223901	BIOPHEN™ UFH Control C2	12 x 1 ml
224001	BIOPHEN™ Arixtra® Control Plasma	2 x 6 x 1 ml
224101	BIOPHEN™ UFH Control C1	12 x 1 ml
224201	BIOPHEN™ LMWH Control C4	12 x 1 ml
224301	BIOPHEN™ LMWH Control Low CI	12 x 1 ml
224401	BIOPHEN™ LMWH Control Low CII	12 x 1 ml
224501	BIOPHEN™ Rivaroxaban Control Plasma	2 x 6 x 1 ml
224701	Dabigatran Control Plasma	2 x 6 x 1 ml
225001	BIOPHEN™ Dabigatran Control Low	2 x 6 x 1 ml
225101	BIOPHEN™ Rivaroxaban Control Low	2 x 6 x 1 ml
225201	BIOPHEN™ Apixaban Control Low	2 x 6 x 1 ml
225301	BIOPHEN™ Apixaban Control	2 x 6 x 1 ml
225401	BIOPHEN™ Edoxaban Control Low	2 x 6 x 1 ml
225501	BIOPHEN™ Edoxaban Control	2 x 6 x 1 ml
226001	BIOPHEN™ Rivaroxaban Calibrator Low	3 x 4 x 1 ml
226101	BIOPHEN™ Apixaban Calibrator Low	3 x 4 x 1 ml
226201	BIOPHEN™ Apixaban Calibrator	3 x 4 x 1 ml
226401	BIOPHEN™ Edoxaban Calibrator Low	3 x 4 x 1 ml
226501	BIOPHEN™ Edoxaban Calibrator	3 x 4 x 1 ml
SC020K	Plasma Hirudin Standard Low	2 x 3 x 1 ml
SC020L	Plasma Hirudin Standard High	2 x 3 x 1 ml
SC025K	Plasma Hirudin control	2 x 3 x 1 ml
SC030K	Argatroban Plasma Calibrator	5 x 4 x 1 ml
SC035K	Argatroban Control Plasma	2 x 6 x 1 ml